

Acrylic Literature Stand

Product Specifications:

Display holds 215mm x 300mm literature in 6 locations.

Foot Print Dimensions: 290mm x 410mm

Product Parts:

- 1 - Acrylic Literature Stand Display Unit
- 2 - Support Hinge
- 3 - Support Hinge Fixture



Instructions:

Step 1 Take all parts out of the case.



Step 2

Lifting from the top side hinges, extend the Acrylic Literature Stand (1) upwards



Step 3

While holding the Acrylic Literature Stand (1) in the upright position, pull the Support Hinges (2) on both sides of the display to face the Support Hinge Fixtures (3). Line up the Support Hinges (2) with the Support Hinge Fixtures (3) and slide the Support Hinge Fixtures (3) over the Support Hinges (2) to secure the Acrylic Literature Stand in its fully extended position.

